


SPE Cartridge Selection Guide


How to Select a Sorbent



One of the most important elements to achieve successful of solid phase extraction is the selection of a sorbent suitable for both the sample matrix and the target analyte.

The sorbent should be carefully selected by taking into account the chemical and physical properties of both the target analyte and the sample matrix. In addition, it is important to develop conditions that are optimal for retaining the target analyte, while removing the sample matrix, then selecting an elution solvent for maximum recovery of the target analyte.

Retentive Capacity of a Sorbent Compared to Sorbent Mass



※Bed volume is the quantity of the solvent necessary to replace the air trapped in the solid phase.
Void volume is equivalent to the bed volume

Recommendation for Selecting an Ion Exchange Sorbent

| Target Analytes | InertSep | pKa* | Structure | Target Ion | | |
|-----------------|-----------------|-----------------------------|------------|--|------------|---|
| | | | | Weak Ion | Strong Ion | |
| Acidic | Anion Exchange | MA-1 4Class Amine | - | -CH ₂ -N ⁺ (R) ₃ | ✓ | × |
| | | MA-2 2Class Amine | 11.0 | -CH ₂ -N (R) ₂ | ✓ | × |
| | | NH ₂ Aminopropyl | 9.8 | -CH ₂ CH ₂ CH ₂ NH ₂ | × | ✓ |
| | | PSA 1Class, 2Class Amine | 10.1, 10.9 | -CH ₂ CH ₂ CH ₂ NHCH ₂ CH ₂ NH ₂ | × | ✓ |
| | | SAX Tri-Methylaminopropyl | - | -CH ₂ CH ₂ CH ₂ N ⁺ (CH ₃) ₃ | ✓ | × |
| | | SAX-2 | - | -CH ₂ CH ₂ CH ₂ N ⁺ (CH ₃) ₃ | ✓ | × |
| Basic | Cation Exchange | MC-1 Sulfonic Acid | 1.0 | -CH ₂ -SO ₃ ⁻ | ✓ | × |
| | | MC-2 Carboxylic Acid | 4.5 | -CH ₂ -COO ⁻ | ✓ | × |
| | | CBA Ethyl Carboxylic Acid | 4.8 | -CH ₂ CH ₂ COO ⁻ | × | ✓ |
| | | PRS Propyl Sulfonic Acid | 1.0 | -CH ₂ CH ₂ CH ₂ SO ₃ ⁻ | ✓ | × |
| | | SCX Benzene Sulfonic Acid | 1.0 | -CH ₂ CH ₂ C ₆ H ₄ SO ₃ ⁻ | ✓ | × |
| | | SCX-2 | - | -CH ₂ CH ₂ C ₆ H ₄ SO ₃ ⁻ | ✓ | × |

* pKa reference value for each functional group.

■ InertSep Series Sorbent Specifications

● Polymer-based Sorbent Specifications

To conduct solid phase extraction, it is necessary to choose the sorbent best suited for the properties of your target compound and sample matrix. The advantages of polymer-based sorbent are the availability in the wide pH range and the absence of secondary interaction which can occur with silica-based sorbents.

| Separation mode | InertSep | Base gel | Functional group | Particle size (µm) | Surface area (m ² /g) | Pore volume (mL/g) | Pore size (nm) | Ion exchange capacity (meq/g) | pH range |
|-----------------|-------------|------------------------|-------------------------------------|--------------------|----------------------------------|--------------------|----------------|-------------------------------|----------|
| Reversed phase | PLS-2 | SDB* ¹ | – | 70 | 700 | 1.1 | 7 | – | 1-14 |
| | PLS-3 | N-MA-SDB* ¹ | – | 60 | 600 | 1.1 | 7 | – | |
| | RP-1 (mini) | MA-DVB* ¹⁺² | – | 70 | 650 | 1.5 | 9 | – | |
| | RP-2 | SDB | weak anion exchanger | 90 | 700 | 0.7 | 4 | – | |
| | Pharma (FF) | N-MA-SDB* ¹ | – | 60 | 600 | 1.1 | 7 | – | |
| Ion exchange | RP-C18 | SDB* ¹ | Octadecyl | 45 | 110 | 0.5 | 18 | – | 1-13 |
| | MA-1 (mini) | MA* ² | Quaternary ammonium | 70 | 250 | 0.7 | 13 | 0.5 | 1-14 |
| | MA-2 (mini) | MA* ² | Diethyl amine | 70 | 250 | 0.8 | 13 | 0.5 | |
| | MC-1 (mini) | MA* ² | Sulfonic acid | 70 | 80 | 0.4 | 20 | 0.5 | |
| | MC-2 (mini) | MA* ² | Carboxylic acid | 70 | 80 | 0.4 | 18 | 0.5 | |
| | MPC | SDB* ¹ | C18, Sulfonic acid | 40 | 100 | – | 18 | – | 1-13 |
| | ME-1 | MA* ² | Iminodiacetic acid | 70 | 80 | 0.5 | 21 | Cu ²⁺ 0.3 mmol/g | 1-14 |
| | ME-2 | MA* ² | Iminodiacetic acid + Tertiary amine | 70 | 80 | 0.5 | 21 | Cu ²⁺ 0.3 mmol/g | |

*1 : In short time, it can be used pH 1 to 14 depending on method.

● Silica-based Sorbent Specifications

The silica-based sorbent materials are more cost-effective and have a higher physical strength compared with polymer-based sorbent materials. Silica offers a wide variety of separation mechanisms using a combination of primary functional group interaction with secondary interactions due to the nature of silica.

| Separation mode | InertSep | Base gel | Functional group | End capped ¹ | Particle size (µm) | Carbon loading (%) | Surface area (m ² /g) | Pore volume (mL/g) | Pore size (nm) | Ion exchange capacity | pH range | | | |
|-----------------|------------------|--------------------------------|----------------------------|-------------------------|--------------------|--------------------|----------------------------------|--------------------|----------------|-----------------------|-------------------|-----|---|-----------------|
| Reversed phase | C18 (FF) | SiO ₂ | Octadecyl (trifunctional) | Excellent | 60 (120) | 19 | 450 | 0.7 | 6 | – | 2-8* ² | | | |
| | C18-B (FF) | | Octadecyl (monofunctional) | Good | 45 (120) | 14 | | | | | | | | |
| | C18-C (FF) | | Octadecyl (trifunctional) | Fair | 60 (120) | 16 | | | | | | | | |
| | C18-ENV | | Octadecyl (trifunctional) | Fair | 60 | 16 | | | | | | | | |
| | C8 | | Octyl | Good | 60 | 12 | | | | | | | | |
| | C8-NE | | Octyl | Poor | 60 | 12 | | | | | | | | |
| | C2 | | Ethyl | Good | 60 | 5.5 | | | | | | | | |
| | CH | | Cyclohexyl | Good | 60 | 7.5 | | | | | | | | |
| | PH | | Phenyl | Good | 60 | 10 | | | | | | | | |
| Ion exchange | SCX | SiO ₂ | Benzenesulfonic acid | None | 45 | 8.5 | 450 | 0.7 | 6 | 0.6 | 2-8* ² | | | |
| | SCX-2 | | Benzenesulfonic acid | None | 60 | 17 | | | | | | 1.2 | | |
| | PRS | | Propylsulfonic acid | None | 45 | 8.5 | | | | | | 1.2 | | |
| | CBA | | Propylcarboxylic acid | None | 45 | 8.5 | | | | | | 1.2 | | |
| | SAX | | Quaternary ammonium | None | 45 | 7 | | | | | | 0.7 | | |
| | SAX-2 | | Quaternary ammonium | None | 60 | 11.5 | | | | | | 0.7 | 6 | 0.45 |
| | PSA | | Ethylenediamine-N- propyl | None | 60 | 11.5 (10.0-13.0) | | | | | | 0.7 | 6 | 1.5 (1.45-1.90) |
| | NH2 | | AminoPropyl | None | 60 | 10 | | | | | | 0.7 | 6 | 0.9 |
| Normal phase | CN | SiO ₂ | Cyanopropyl | None | 45 | 0.7 | 450 | 0.7 | 6 | – | 2-8* ² | | | |
| | 2OH | | Diol | None | 60 | 10 | | | | | | | | |
| | Si | | – | None | 60 | – | | | | | | | | |
| | AL | Al ₂ O ₃ | Aluminium oxide | None | 100 | – | 130 | 0.3 | 8 | – | | | | |
| | FL | MgO·SiO ₂ | Magnesium silicate | None | 50-200 | – | 230 | 0.5 | 9 | – | | | | |
| FL-PR | SiO ₂ | | None | 100-300 | – | 230 | 0.5 | 9 | – | | | | | |

*1 : Styrene divinylbenzene copolymer

*2 : Methacrylate polymer

● Specialty Phases

| InertSep | Base gel | Particle size | Surface area (m ² /g) | Pore volume (mL/g) | Pore size (nm) |
|----------|------------------|---------------|----------------------------------|--------------------|----------------|
| GC | Graphite Carbon | 120/400 mesh | 85 | 1 | 45 |
| GC-e | Graphite Carbon | 100/200 mesh | 90 | 1 | 50 |
| AC | Activated Carbon | 65/150 mesh | 800-1200 | – | – |