



# MonoSpin

Solid Phase Extraction Spin Column




## Contents

About MonoSpin .....	2
Product line-up .....	4
Features of MonoSpin .....	5
Application .....	7
Spin column for Phospholipid	
MonoSpin Phospholipid .....	13
Spin column for purification of antibody	
MonoSpin ProA , ProG .....	16
MonoSpin 96-well plate .....	20
Accessory .....	21
Ordering Information .....	22

# Silica monolith ~ New separation media that are neither particulate nor membrane~

Silica monoliths are integral silica gels with uniform continuous pores synthesized from ethyl silicate. Unlike the particle media, Silica monolith is shaped like a disk. Silica monoliths have high liquid permeability and large surface area as it has structure having through-pores and mesopores on the framework surface. Therefore, this state-of-the-art media is becoming popular globally for its characteristics: high recovery, high performance of adsorption and desorption.




## Different points from the particle packing type pretreatment columns

Disk-shaped silica monoliths isn't required to use frits to hold particle media like conventional SPE cartridges. It also has a very large surface area, which makes it possible to reduce the volume of sample. From this point, Silica monolith makes it possible to prevent from remaining samples in the cartridge as well as complete elution for small amount of sample.

Despite the fact that it has high liquid permeability, It is also suitable for rapid elution without losing its high recovery as it enables fast sample diffusion and quick separation.


### Silica monolith



- No filter required
- Minimized separation media

Bed volume for separation media : **small**  
Sample diffusion in the column : **fast**  
Separation Speed : **fast**

### Particle-filled Form




- Need for filters  
→ liquid may be remained in the filter


Bed volume for separation media : **large**  
Sample diffusion in the column : **slow**  
Separation Speed : **slow**

# MonoSpin

MonoSpin is a SPE spin-column using silica monoliths with uniform continuum pores. It effectively and quickly extracts, isolates, purifies and concentrates samples by centrifugation with an ease.



MonoSpin




Silica monolith

## Feature

- Easy operation by centrifuge
- Speedy sample treatment with a superb through pore
- Excellent reproducibility (S-type) even at elution volumes of 100  $\mu$ L or less.

## Operation method

Whole process for sample treatment can be done within 10 minutes.




## Shape

MonoSpin series cartridge with different types are available:

Type S: Excellent for the pretreatment of the sample for 50-800 $\mu$ L

Type L: Suitable for sample 0.5-8mL.

For the details of the varied functional group, please see the next page.



### S Type

- Disk size:  $\phi 4.2 \times 1.5$  mm
- Sample volume up to 800  $\mu$ L
- Elution volume: 50 to 800  $\mu$ L
- Centrifugation speed: 2,000 to 10,000 x g




### L Type

- Disk size:  $\phi 9 \times 3$  mm
- Sample volume to 8 mL:
- Elution volume : 0.5-8 mL
- Centrifugation speed : 1,000 x g

(NOTE) The MonoSpin ProA, MonoSpin ProG has different shapes. Please see page 16 for details.

### 96 Well plate type




- Sample volume up to 800  $\mu$ L
- Elution volume :50 ~ 800  $\mu$ L
- Centrifugation speed : 1,000 to 5,000  $\times$  g (can be used in vacuum aspiration)

(NOTE) The MonoSpin C18 FF, MonoSpin ProA, MonoSpin ProG has different specifications. Please see page 14, 15 for details.

# MonoSpin Series Lineup

## MonoSpin C18/C18 FF

**S** **L** **96**




\* Only S-type for MonoSpin C18 FF.

Octadecyl functional group. Optimal for drug extraction in biological samples, and desalting & enrichment of peptide samples. High-flow (FF) specification is also available.

## MonoSpin C18-AX


**S** **96**



Is a mix mode type in which both octadecyl and quaternary ammonium groups are chemically bonded. It can reliably retain bio samples in high salt concentration. It is especially suitable for recovery of acidic drugs.

## MonoSpin SAX


**S** **L** **96**



Bonded with Trimethyl aminopropyl combining both strong anion exchange & weak hydrophobic interaction. Optimal for the extraction of acidic drugs.

## MonoSpin NH2


**S** **L** **96**



Bonded with aminopropyl. Optimal for the enrichment of sugar chain and/or hydrophilic compounds by HILIC mode.

## MonoSpin Amide

**S** **96**



Bonded with amide group. Optimal for the extraction of sugar chains and various acidic and basic hydrophilic compounds by HILIC mode.

## MonoSpin TiO


**S**



Monolith skeleton coated with dioxide titanium. Excellent for the enrichment of phosphopeptides

## MonoSpin ME


**S** **L** \*Please refer to P12



Bonded with iminodiacetic acid groups. It is optimal for the recovery of trace metals in samples.

## MonoSpin ProA


**S** **96** \*Please refer to P15



Protein A is immobilized on the Monolith. It enables you efficient purification of Antibodies.

## MonoSpin Ph


**S**



Phenyl group is chemically bonded which makes it possible to utilize weaker hydrophobicity than C18. It is suitable for the recovery of hydrophobic drugs from biological samples under reversed phase mode.

## MonoSpin C18-CX


**S** **96**



Octadecyl and Benzenesulfonic acid groups are bonded. It is suitable for the purification of dissociated basic drugs in serum and urine. Compared with MonoSpin C18 and SCX alone, SCX has higher cleanup efficacy as it works as Hydrophobic and ion-exchange interactions.

## MonoSpin SCX


**S** **L** **96**



Bonded with propyl benzenesulfonic acid combining both strong cation exchange & hydrophobic interaction. Optimal for the extraction of basic drugs.

## MonoSpin CBA


**S** **L** **96**



Bonded with propyl benzenesulfonic acid combining both strong cation exchange & hydrophobic interaction. Optimal for the extraction of basic drugs.

## MonoSpin PBA


**S** **96**



Bonded with phenyl boric acid which gives you higher selectivity. Excellent for the extraction of cis diol compounds, such as catechol amines.

## MonoSpin Trypsin

**S**



Columns immobilized with trypsin, a protein digestive enzyme. It enables rapid digestion of proteins.

## MonoSpin Phospholipid


**S** **L** \*Please refer to P13



Phospholipid Removal Column coated with titanium dioxide and zirconium dioxide on Silica monolith. It adsorbs phospholipids in samples with an easy pretreatment.

## MonoSpin ProG

**S** **96** \*Please refer to P15



Protein G is immobilized on the Monolith. It enables you efficient purification of Antibodies.

**S** :S-type column products.

**L** :L-type column products.


**96** :96-well plate-type product.

# Characteristics of MonoSpin Series


## Purification and enrichment of trace analytes

Due to its high permeability, MonoSpin series enable you more faster and efficient purification and enrichment with centrifugation. It's also recommended for the elution of small volume samples, and trace analytes can be collected without dilution.

### S Type



### L Type



## Physical properties of MonoSpin series

Product name	Functional Group	S/96 well		L Type		Surface Area (m <sup>2</sup> /g)	Bed Capacity (for type S).
		Through pore (μm)	Mesopore (nm)	Through pore (μm)	Mesopore (nm)		
MonoSpin C18	Octadecyl group	5	10	10	10	350	100 μg (amitriptyline)
MonoSpin C18 FF	Octadecyl group	20	15	10	10	300	50 μg (amitriptyline)
MonoSpin Ph	Phenyl group	5	10	-	-	350	100 μg (amitriptyline)
MonoSpin C18-AX	Octadecyl group, Quaternary ammonium	5	10	-	-	350	100 μg (ibuprofen)
MonoSpin C18-CX	Octadecyl group, Benzenesulfonic acid group	5	10	-	-	350	100 μg (amitriptyline)
MonoSpin SAX	Trimethylaminopropyl group	5	10	10	10	350	100 μg (ibuprofen)
MonoSpin SCX	Propylbenzenesulfonic acid group	5	10	10	10	350	100 μg (amitriptyline)
MonoSpin NH2	Aminopropyl-group	5	10	10	10	350	100 μg (maltopentaose)
MonoSpin CBA	Carboxyl group	5	10	10	10	350	100 μg (amitriptyline)
MonoSpin Amide	Amide group	5	10	-	-	350	100 μg (angiotensin II)
MonoSpin PBA	Phenylboronic acids	5	10	-	-	350	100 μg (dopamine)
MonoSpin TiO	Titanium dioxide	20	15	-	-	350	40 μg (adenosine monophosphate)
MonoSpin Trypsin	Trypsin	5	10	-	-	350	-
MonoSpin ME	Iminodiacetic acid group	5	10	10	10	350	25 μg (Cu ions)
MonoSpin Phospholipid	Titanium dioxide Zirconium dioxide	5	10	10	10	350	10 μL (human serum)
MonoSpin ProA	Protein A	2	60	-	-	-	400 μg (human IgG)
MonoSpin ProG	Protein G	2	60	-	-	-	400 μg (human IgG)

## Specifications for Shape and Type

Type	MonoSpin S type <sup>*1</sup>	MonoSpin FF <sup>*2</sup>	MonoSpin L type	MonoSpin 96 well type
Disk size	Φ4.2 × 1.5 mm	Φ4.2 × 1.5 mm	Φ9 × 3 mm	Φ4.2 × 1.5 mm
Sample Volume	Up to 800 μL	Up to 800 μL	Up to 8 mL	Up to 800 μL
Elution Volume	50~800 μL	50~800 μL	0.5 ~ 8 mL	100~800 μL
Centrifugal force	2,000~10,000 × g	1,000 × g	1,000 × g	1,000~5,000 × g


\* 1: MonoSpin ProA and MonoSpin ProG are different in specifications. Please refer to page 15 for the details.

\* 2: FF type is available for MonoSpin C18 FF only.

## Which MonoSpin is suitable for your analyte?

The MonoSpin series are optimized as a spin-column for pretreatment of biological samples. If you are working on highly viscous samples such as blood, MonoSpin C18 FF is the best choice.

Please refer to the following chart for choosing appropriate MonoSpin.






# Application

## Purification of Amphetamines in urine using MonoSpin C18

Sample Volume: 800  $\mu$ L

Urine: 400  $\mu$ L  
 Buffer solution (pH 13): 400  $\mu$ L  
 \* Sample was mixed for 1 min at 10,000 x g. Transferred and used the supernatant as sample.



### 1. Conditioning

- ① Methanol 300  $\mu$ L
- ② Buffer (pH 13) 300  $\mu$ L

### 2. Adsorption

Sample solution 800  $\mu$ L

### 3. Wash


Buffer (pH 13). 300  $\mu$ L

### 4. Elution

Methanol-0.1% Formic acid (1:1,v/v) 100  $\mu$ L

Purified sample

Centrifuge : 5,000  $\times$  g



### Conditions


Column : InertSustainSwift C18 (3  $\mu$ m, 150  $\times$  2.1 mm I.D.)  
 Eluent A) 10 mM ammonium acetate-formic acid (pH 3.3).  
 B) CH<sub>3</sub>OH  
 A/B = 90/10 - 2 min - 90/10 - 13 min - 70/30,v/v  
 Flow Rate : 0.3 mL/min  
 Col. Temp. : 40  $^{\circ}$ C  
 Detection : LC/MS  
 Sample : 1. Norephedrine  
 2. Ephedrine  
 3. Methylephedrine  
 4. Amphetamine  
 5. Methamphetamine  
 6. 3,4-methylenedioxyamphetamine  
 7. 3,4-methylenedioxymethamphetamine

※ Data provided by Dr. Namera, Hiroshima University

## Recovery of drugs in biological samples using MonoSpin C18<sup>®</sup>

Sample Volume 600  $\mu$ L

Serum: 200  $\mu$ L  
 10 mM potassium phosphate:  
 400  $\mu$ L (pH 7.0).  
 \* Sample was mixed for 1 min at 10,000 x g. Transferred and used the supernatant as sample.



### 1. Conditioning

- ① Methanol 300  $\mu$ L
- ② 10 mM potassium phosphate, 300  $\mu$ L (pH 7.0).

### 2. Adsorption

Sample solution 600  $\mu$ L

### 3. Wash

Water 300  $\mu$ L

### 4. Elution

Acetonitrile 200  $\mu$ L

Purified sample

Centrifuge : 2,300  $\times$  g

Day-to-day reproducibility of drug in serum using MonoSpin C18 (3 days, n = 10).

Sample	Concentration (ng/mL)	Recovery rate (%)	RSD (%)
Desipramine	5	91.2	4.8
	10	86.1	3.3
	50	85.2	5.9
	250	88.4	6.5
Imipramine	5	96.3	9.5
	10	95.8	1.5
	50	94.5	0.9
	250	95.9	0.9
Fluvoxamine	5	96.8	11.6
	10	87.1	5.0
	50	86.8	8.1
	250	87.5	9.7


Sample	Concentration (ng/mL)	Recovery rate (%)	RSD (%)
Paroxetine	5	83.7	3.9
	10	84.1	7.8
	50	83.9	8.2
	250	86.7	7.5
Maprotiline	5	85.7	8.1
	10	84.7	3.2
	50	88.6	5.4
	250	87.5	7.7
Duloxetine	5	106.3	9.9
	10	104.8	6.7
	50	99.8	8.7
	250	99.8	6.0

Sample	Concentration (ng/mL)	Recovery rate (%)	RSD (%)
Amitriptyline	5	83.7	7.0
	10	81.8	2.8
	50	83.8	3.0
	250	88.4	2.7
Sulpiride	5	97.9	9.0
	10	95.5	8.5
	50	90.8	2.6
	250	92.6	3.0

## Desalination of protein digestion using MonoSpin C18<sup>®</sup>

Maximum sample solution=800 µL

After Tryptic digestion, add TFA to adjust the concentration to 0.1%.



Centrifuge : 2,300 × g

### 1. Conditioning

Acetonitrile 200 µL  
→ 0.1% aqueous TFA  
200 µL

### 2. Adsorption

Sample solution  
Maximum 800 µL


### 3. Wash

0.1% aqueous TFA  
200 µL


### 4. Elution

60% acetonitrile  
200 µL


### Desalted samples




Before purification



After purification



Enlarged view



Enlarged view

### Conditions

Column :Inertsil ODS-3  
(3 µm, 150 × 2.1 mm I.D.)  
Eluent :A)H<sub>2</sub>O (0.1 % TFA)  
B) Acetonitrile (0.1 % TFA)  
A/B = 90/10 - 20 min - 50/50  
Flow Rate :UV 210 nm  
Col. Temp. :0.2 mL/min  
Detection :40 °C  
Sample :Digested BSA 2 µL

Highly concentrated denaturant and salt in digestive were successfully removed using MonoSpin C18.

## Rapid Digestion of BSAs by MonoSpin's Trypsin HP


### ■ Ex. Reductive alkylation protocol

1 mg bovine serum-albumin

- 500 mM Tris-HCL (pH 8. 0)-- 8M urea (Solution 1): 175µL
- 40 mg/mL Dithiothreitol in Solution 1: 25µL
- Incubation at 37 °C for 90 min
- 40 mg/mL Iodoacetamide in Solution 1: 50µL
- Incubation at 37 °C for 30 min (under shaded conditions)

Reductive alkylation of proteins: 250µL

- Dilute with 50mM Ammonium bicarbonate to adjust the urea to 2M: 750µL




### Conditions

Column :Inertsil ODS-3  
(3 µm, 150 × 2.1 mm I.D.)  
Eluent :A)H<sub>2</sub>O (0.1 % HCOOH)  
B) Acetonitrile (0.1 % HCOOH)  
A/B = 90/10 - 20 min - 50/50  
Flow Rate :UV 210 nm  
Col. Temp. :0.2 mL/min  
Detection :40 °C  
Sample :Digested BSA 2 µL


### MonoSpin Trypsin HP

NOTE) The method of reductive alkylation should be optimized depending on the type of protein.

#### ● In-Solution digestion (37°C for 10 h)



#### ● Digested with MonoSpin Trypsin HP (at 25°C for 10 min)



Trypsin-immobilized spin column can complete the process just in 10 min.

NOTE) For digestion, be sure to use protein after reductive alkylation.




# Application

## Fractionation of Protein digests using MonoSpin SCX.

The use of spin columns and elution salt concentration in stepwise makes it possible to fractionate peptides without using 2D-LC systems or other complex systems.

Sample Volume: 500  $\mu$ L

Used Peptide sample dissolved in 0.1% Formic acid after desalting with MonoSpin C18.



Apply the eluent, centrifuge, and then attach a new tube to apply the next eluent.


Each eluate composition

- ① 25 mM HCOONH<sub>4</sub> 200  $\mu$ L
  - ② 50 mM HCOONH<sub>4</sub> 200  $\mu$ L
  - ③ 100 mM HCOONH<sub>4</sub> 200  $\mu$ L
  - ④ 500 mM HCOONH<sub>4</sub> 200  $\mu$ L
  - ⑤ 1 M HCOONH<sub>4</sub> 200  $\mu$ L
- Injection) Each solution contains 10% acetonitrile.

### Conditions

Column :Inertsil ODS-3 (3  $\mu$ m, 2.1  $\times$  150 mm)  
Eluent :A)H<sub>2</sub>O (0.1% HCOOH)  
B) Acetonitrile (0.1% HCOOH)  
A/B = 90/10 - 20 min - 50/50

Detection :UV 210 nm  
Flow Rate :0.2 mL/min  
Col. Temp. :40  $^{\circ}$ C  
Injection Vol. :2  $\mu$ L




It is possible to perform the fractionation of peptide samples by applying eluent with different salt concentrations.

## Purification of pyridylaminated glycans using MonoSpin's NH<sub>2</sub>.


Sample volume: 800  $\mu$ L

Dissolve the sample to adjust the concentration of ACN to 90~95%.




Centrifuge 2,300  $\times$  g

Before purification



Purified with MonoSpin NH<sub>2</sub>



### Conditions


Column :NH<sub>2</sub> Column (5  $\mu$ m, 250  $\times$  4.6 mm I.D.)  
Eluent :A)H<sub>2</sub>O/Acetonitrile  
=5/95 0.1% Formic acid  
B) H<sub>2</sub>O/Acetonitrile  
=95/5 0.1% Formic acid  
A/B = 90/10-10 min-90/10-40 min-60/40  
Flow Rate :1 mL/min  
Detection :FL Em 320 nm, Ex 400 nm  
Injection Vol. :1.5  $\mu$ L

## Purification of Paraquat and Diquat using MonoSpin CBA®

Sample volume 600 µL

Urine: 200 µL  
10 mM potassium phosphate  
Buffer solution (pH 7.0): 400 µL




Centrifuge : 10,000 × g



Standard Solution (1 µg/mL)

Urine + pesticide (1 µg/mL each)

After purification with MonoSpin CBA



### Conditions



Column : Inertsil ODS-3 (5 µm, 150 mm × 4.6 mm I.D.)  
Eluent : 0.2 M phosphoric acid, 0.1 M diethyl amine, 7.5 mM IPCC08(IPCC-0.8, Sodium 1-Octanesulfonate) /Acetonitrile=89/11  
Flow Rate : 1 mL/min  
Col. Temp. : 40 °C  
Detection : PDA 290 nm  
Injection Vol. : 50 µL

## Recovery of hormones in serum using MonoSpin C18

Sample preparation

Add 10 µL of 1 mg/mL adrenomedullin to serum: 190 µL.  
Centrifuge the sample after addition of 0.1% TFA solution 200 µL.  
Used the supernatant as sample.

Centrifuge : 2,300 × g



### Conditions

Column : InertSustain C18 (2 µm, 50 × 2.1 mm I.D.)  
Eluent : A) 0.1 % TFA in Water  
B) 0.1 % TFA in Acetonitrile  
A/B = 85/15 – 5 min – 50/50 – 2 min-50/50  
Flow Rate : 200 µL/min  
Col. Temp. : 40 °C  
Detection : UV 210 nm  
Injection Vol. : 10 µL


# Application

## Purification of Catecholamines using MonoSpin PBA.

Sample solution 250  $\mu$ L




Sample solution (urine or serum)  
200  $\mu$ L  
1 M Dipotassium hydrogen phosphate  
(adjusted to pH 8 with phosphoric acid)  
50  $\mu$ L

Centrifuge : 10,000  $\times$  g



Before Purification

Purification with MonoSpin PBA (sample amount: 500  $\mu$ L)



### Conditions

Column : Inertsil ODS-3  
(5  $\mu$ m, 150 mm  $\times$  2.1 mm I.D.)  
Eluent : 50 mM Phosphate Buffer (pH 5.6)  
50 mg/L EDTA  
600 mg/L IPCC-008  
-10 % Methanol  
Flow Rate : 0.3 mL/min  
Col.Temp. : 35  $^{\circ}$ C  
Injection : 5  $\mu$ L  
Detection : ECD Pulse Mode  
Sample : 1. Noradrenaline  
2. Adrenaline  
3. DHBA  
4. Dopamine


By using MonoSpin PBA, we can selectively recover and purify compounds with cis-type diols such as catecholamines. See our website Technical Note LT093 for more information.

## Purification of Organophosphorus pesticides in human serum using MonoSpin TiO

Sample Volume 50  $\mu$ L



Sample 10  $\mu$ L  
Water 40  $\mu$ L

Centrifuge : 5,200  $\times$  g



Before purification






After purification using MonoSpin TiO



Bialaphos  
Glyphosate  
MPPA  
Glufosinate  
AMPA

### Conditions

Column : ODS Column (150  $\times$  2.1 mm I.D.)  
Eluent : A) CH<sub>3</sub>OH  
B) 20 mM HCO<sub>2</sub>NH<sub>4</sub> (pH 3.0)  
A/B = 15/85, v/v  
Flow Rate : 200  $\mu$ L/min  
Detection : SIM  
Injection Vol. : 5  $\mu$ L  
Sample : 1. Bialaphos  
2. Glyphosate  
3. MPPA  
4. Glufosinate  
5. AMPA






MonoSpin TiO shows selectivity for phosphate sites in compounds.

# Application

## With Insulin or BSAs

### Sample Preparation

Adjust concentration of Insulin and BSA with 0.1% aqueous TFA.





Please see Technical Note LT157 for more information .

## Analysis of blood sample using MonoSpin C18FF

### Sample Preparation

Mix blood sample (0.3mL) and 300mM phosphate buffer (pH 10.0).

Use supernatant as sample after centrifugation at 12,000 x g for 5 min.



### Conditions

Column: InertSustain Phenyl (3  $\mu$ m, 150  $\times$  2.1 mm I.D.)  
 Eluent: acetonitrile-HCOONH<sub>4</sub> (10 mM, 0.1% HCOOH) = 25:75 (v/v)  
 Flow Rate: 0.2 mL/min  
 Col. Temp.: 40  $^{\circ}$ C  
 Detection: MS(ESI)

# MonoSpin Phospholipid

Monolith



ZrO<sub>2</sub>


Phospholipid Removal Column coated with titanium dioxide and zirconium dioxide on Silica monolith. It adsorbs phospholipids in samples such as blood and serum with an easy pretreatment. More importantly, the adsorbed phospholipids can also be recovered very well.

Cartridge shape: S-type, L-type

Functional groups: titanium dioxide, zirconium dioxide


## Features

- Phospholipids in biological sample can be removed with an easy operation.
- It decreases matrix effect significantly since it removes more than 90% of Phospholipids.
- Perform well even for small volume sample
- The adsorbed phospholipids can also be recovered.



## How it's adsorbed?


- Monolith skeletal structure coated with TiO<sub>2</sub> and ZrO<sub>2</sub> Selectively interacts with metal oxides and phosphorylated compounds, resulting in removing more than 90% of phospholipids!



## Operation

### Sample Preparation

Mix 0.1 % formic acid and serum sample 4:1 into 2mL micro tube.  
 ↓  
 Use supernatant after centrifugation at 10,000 × g 30 sec



### Related Product: FastRemover® for Phospholipid



The FastRemover® for Phospholipid 96-well plate deliver a rapid and effective removal of proteins and phospholipids in plasma and serum samples without sacrificing the recovery of your target analytes.

# Publicly Available Reference

Functional Group	Compounds	Reference
C18	Pyrrolidinophenone type designer	J. Chromatogr. B, 2013, 30, 942-943
	Aconitines and Colchicine	Chromatographia, 2015, 78(15), 1041-1048
	Eperisone, Tolperisone, and Tizanidine	J.AOAC Int. 2014, 97(6), 1546-1551
	MAM-2201	Forensic Toxicology, 2013, 31(2), 333-337
	Diquat, Paraquat	Anal. Bioanal. Chem., 2011, 400(1), 25-31
	Nanoparticles	J. Chromatogr. A, 2015, 1404, 141-145
	Iodide	Am. J. Mod. Chromatogr., 2015, 2(1), 1-6
	$\alpha$ -Pyrrolidinovalerophenone	Forensic Toxicol., 2014, 32(1), 68-74
	Organophosphorus compounds	Anal. Sci., 2011, 27(10), 999-1005
	review	Bioanalysis., 2015, 7(17), 2171-2176
	Phthalate esters	J. Pharm. Anal., 2011, 1(2), 92-99
	Oxidized phospholipids	J. Lipid. Res., 2017, 58(11), 2229-2237
	Medicinal toxicants	J. Clin. Pharm. Ther., 2017, 42(4), 454-460
	N-1-Naphthalenyl-1-pentyl-1H-indole-3-carboxamide	Forensic Toxicol., 2015, 33(1), 165-169
	Peptides	Cancer Res., 2017, 77(4), 926-936
	Peptides	Bio protocol, 2015, 5(8), 2015
	oxPUFAs	Sci. Rep., 2018, 8, 7954
	Desalting	Amino Acids., 2018, 50(1), 117-124
	Peptides	Clin. Exp. Nephrol., 2018, 22(4), 782-788
	25-Hydroxyvitamin D3	Anal. Sci., 2018, 34(9), 1043-1047
	Peptides	Biosci. Biotechnol. Biochem., 2017, 81(12), 2237-2243
	Desalting	Org. Biomol. Chem., 2018, 17(1), 165-171
	Peptides	Methods Mol. Biol. 2018, 1696, 91-105
	Purines	Nucleosides Nucleotides Nucleic Acids, 2018, 37(6), 348-352
	Flavonoid	J. Chem. Ecol. 2016, 42(12), 1226-1236
	Peptides	Biosci. Biotechnol. Biochem., 2018, 82(8), 1309-1315
	Peptides	Data Brief., 2018, 31(17), 604-609
	Desalting	J. Proteomics, 2018, 181, 238-248
	Peptides	Data Brief., 2017, 12(11), 252-257
	glucocorticoids	J. Chromatogr. B, 2017, 1057, 62-69
	Peptides	Bioresour. Technol., 2018, 254, 278-283
	ITRAQ labeled desalting	Int. J. Oncol., 2015, 47(1), 384-390
	Peptides	Biomass Bioenergy, 2016, 91, 83-90
	Peptides	Neurogenetics, 2019, 20(1), 9-25
	Desalting	J. Pept. Sci., 2018, 24(12), e3133
	Desalting of LaIT1	Mass Spectrometry, 2017, 6(1), A0059
	Peptides	J. Proteomics, 2015, 119, 183-195
	Liraglutide	J. Chromatogr. B, 2018, 109, 29-35
	Desalting of LaIT1	J. Pept. Sci., 2015, 21(8), 636-643
	Peptides	Proc. Natl. Acad. Sci., 2018, 115(14), 3646-3651
	Peptides	Oncogene, 2017, 36(26), 3740-3748. doi: 10.1038/onc.2016.524
	Plant samples	Sci. Rep., 2017, 7(1), 1243. doi: 10.1038/s41598-017-01390-3
	Peptides	Sci. Rep., 2018, 22, 8(1), 1303
	Peptides	Sci. Rep., 2016, 6, 26723
	Drugs	J. Chromatogr. B, 2008, 867(1), 99-104
	Dibucaine, Naphazoline	J. Chromatogr. B, 2008, 872, (1-2), 186-190
	Amphetamines	J. Chromatogr. A, 2008, 1208(1-2), 71-75
	Drugs	Chromatographia, 2009, 70(3), 519-526
	Amphetamines	Anal. Chim. Acta, 2010, 661(1), 42-46
	Eperisone, Tolperisone	J.Health Sci., 2010, 56(5), 598-605
Diquat, Paraquat	Anal. Bioanal. Chem., 2011, 400(1), 25-31	
[1-(5-fluoropentyl)-1H-indol-3-yl](4-methyl-1-naphthalenyl)methanone (MAM-2201)	Forensic Toxicol., 2013, 31(2), 333-337	
a-Pyrrolidinovalerophenone (a-PVP) and a-pyrrolidinobutiophenone (a-PBP)	Forensic Toxicol., 2014, 32, 68-74	
Pyrrolidinophenone-type designer drugs	J. Chromatogr. B, 2013, 942-943, 15-20	
Phthalates	J. Pharm. Anal., 2011, 1(2), 92-99	
Peptides	Proteomics, 2013, 13(5), 751-755	
Peptides	J. Proteomics., 2013, 84(12), 40-51	
Naringin	J. Clin. Pharmacol., 2013, 53(7), 738-745	



# Publicly Available Reference

Functional Group	Compounds	Reference	
C18 FF	Drugs	J. Chromatogr. A, 2017, 1517, 9-17	
C18, C18CX	Cardiovascular drug	Acta Chromatographica, <a href="https://doi.org/10.1556/1326.2018.00493">https://doi.org/10.1556/1326.2018.00493</a>	
C18, SCX	Melamine	J. Anal. Sci. Meth. Instrum., 2012, 2, 68-73	
C18, SCX	Peptides	Sci. Rep., 2017, 7(1), 11137	
C18, TiO	Peptides	Int. J. Mol. Sci., 2018, 19(9), 2655	
C18,SAX	Aamphetamines, Opiates, and THC	Forensic Toxicol., 2013, 31(2), 312-321	
C18AX	Oxidized Fatty Acids	Mod. Chem. Appl., 2015, 3, 3	
C18-CX	Clean up	J. Occup. Health., 2018, 60(2), 140-147	
	Arsenite, Arsenate, and Methylarsenate	J. Sep. Sci., 2012, 35(18), 2506-2513	
	Drugs	J. Sep. Sci., 2011, 34(16-17), 2232-2239	
	Halogenated compounds	Toxicology, 2013, 314(1), 22-9	
CBA	clenbuterol	Talanta, 2018, 186, 521-526	
CBA, Amide	Tetrodotoxin	Chromatographia, 2014, 77, (9-10), 687-693	
Amide	PA-labelled glycans	Bicsci. Biotechnol. Biochem., 2012, 76(10), 1982-1983	
NH2	PA labeled N-glycans	Glycoconj. J., 2017, 34(4), 537-544	
	PA-labelled glycans	Plant Biotechnol. J., 2016, 14(8), 1682-1694	
	Oligosaccharides	Sci Rep. 2017, 26(7) :46099. doi	
	Pyridylaminated Oligosaccharides nanoparticles	Anal. Sci., 2016, 32(5), 487-490 J. Sep. Sci., 2015, 38, 283-290	
PBA	Pyridylamino monosaccharide	Bicsci. Biotechnol. Biochem., 2011, 75(7), 1405-1407	
	Catecholamines	J. Comp. Neurol., 2016, 524(18), 3849-3864	
	Catecholamines	Food Chem., 2019, 276, 376-382	
	Catecholamines	EBioMedicine., 2016, 8, 60-71	
	Catecholamines	PLoS One, 2018, 13(7), e0201203	
	hippocampal monoamines	J. Pharmacol. Sci., 2016, 132(4), 249-254	
	Catecholamines	J. Chromatogr. B, 2015, 985, 142-148	
	Catecholamines	Biol. Pharm. Bull., 2017, 40(2), 227-233	
	Catecholamines	Biosci. Biotechnol. Biochem., 2018, 82(3), 497-506	
	Serotonine and Noradrenaline	Br. J. Pharmacol., 2015, 172(5), 1250-1262	
Phospholipid	Cis-diol groups	Anal. Chim. Acta., 2015, 857(1), 64-70	
	Allergenic ingredients	Food Control, 2018, 84, 89-96	
	Adenosine	Biosens. Bioelectron., 2013, 15(41), 379-385	
	Farnesyl pyrophosphate	Anal. Bioanal. Chem., 2017, 409(14), 3551-3560	
	ProteinA, G	IgG	Biochimie., 2018, 145, 113-124
ProteinG	IgG	Virology, 2019, 15, 527, 132-140	
	IgG	PLoS One, 2017, 12(7):e0181181	
	IgG	Bioanalysis, 2018, 10(18), 1501-1510	
SAX	Deoxyribonucleoside	Biotechnol., 2016, 228, 52-57.	
	Alendronate	Legal. Medicine, 2018, 30, 14-20	
	Urinary excretion	Nucleosides Nucleotides Nucleic Acids. 2016, 35(10-12), 559-565.	
	metabolite of 18 F-THK5351	Eur. J. Nucl. Med. Mol. Imaging, 2016, 43(12), 2211-2218	
SCX	Methylated lysine	Anal. Bioanal. Chem., 2018, 410(17), 4189-4194	
	Amino acid	Psychiatry Res., 2016, 238, 203-210	
	Amino acid	J. Sep. Sci., 2014, 37(16), 2087-2094	
	Angiogenic peptide	BioSci. Trends, 2016, 10(6), 500-506	
	iTRAQ-labeled peptides	Biochim. Biophys. Acta, 2018, 1865(6), 874-888	
	Amino acid	Sci. Rep., 2018, 8(1), 14587	
	Morphine, Codeine, Dihydrocodeine	J. AOAC Int., 2011, 94(3), 765-774.	
TiO	Fluorescence derivatization	Biomed. Chromatogr., 2012, 26(2), 147-151	
	Amino acid	Orig. Life Evol. Biosph., 2013, 43(2), 99-108	
	Glyphosate	Acta Chromatographica, <a href="https://doi.org/10.1556/1326.2018.00513">https://doi.org/10.1556/1326.2018.00513</a>	
	Trypsin	Protein digestion	J. Am. Chem. Soc., 2018, 140(38), 11982-11991
		Protein digestion	Anal. Sci., 2018, 34(4), 397-406
review		Forensic Toxicol., 2010, 28(2), 61-68	
review		Trac. Trends Anal. Chem., 2013, 45, 182-196	
review		Electrophoresis. 2017, 38(22-23), 2851-2869	
review	review	Chromatogr., 2015, 2(1), 79-95	
	review	J. Pharm. Biomed. Anal., 2018, 161, 51-60	


# Rapid Purification of Antibodies(1)- MonoSpin ProA, MonoSpin ProG

MonoSpin® ProA and MonoSpin® ProG are available already immobilized onto a silica monolith offering rapid purification of antibodies. A 96-well plate format is available to purify a multi-analyte. Each reagent for purification of samples is attached.



## Features

The silica is modified with a hydrophilic polymer and then immobilized with either Protein A or Protein G to prevent the adsorption of proteins, resulting in higher purification and recovery of antibodies.



Silica monolith surfaces immobilized with Protein A, Protein G have modified hydrophilic polymers, which suppress the nonspecific adsorption of proteins and allow the recovery of more pure antibodies.

Specification	
Bonded phase	Protein A or Protein G
Through-pore size	2 μm
Meso-pore size	60 nm
Column size	Φ 4.2 × 1.5 mm
Sample Volume	50-500 μL
Centrifugation speed	2,300 ×g *
Recovery rates	MonoSpin ProA IgG 90%(With 400μg IgG)
	MonoSpin ProG IgG 90%(With 300μg IgG)

\* :96-well plate type can also be used with vacuum aspiration (e.g. -0.015 MPa).

Silica Monolith is available for different shapes

### Spin Column Type




- Purification with compact tabletop centrifuge just in two minutes(e.g. 2,300 × g)
- Suitable for purification of small volume sample(up to 0.4mg)

### 96 Well plate type



- Purification by both aspiration or centrifuge
- Available for a multi-analyte with same spin column volume..

### Large Spin Column




- Max. 16mg antibody can be recovered by centrifuge.



# Rapid Purification of Antibodies(2)- MonoSpin ProA, MonoSpin ProG

## Purification of IgG Using MonoSpin® ProA and MonoSpin® ProG in Only 5 min.

Sample can be neutralized and maintained in stable condition by putting neutral solution in tube when recovering anti body.




As shown below, the antibody concentrations were determined quantitatively from medium of CHO cells. The purified antibodies show very less impurities by the results from electrophoresis.



## Enrichment of Antibody Solution Using MonoSpin® ProA


500  $\mu$ L volume of 0.025 mg / mL of human IgG solution is applied to MonoSpin® ProG spin column for 10times (In = I1 – I10). And then the elution of IgG concentration is measured with 100  $\mu$ L elution buffer twice (En = E1 and E2). The first IgG elution (E1) is 50 fold concentration of standard solution and indicates 90 % recovery of IgG without the loss of IgG.




# MonoSpin ProA, MonoSpin ProG

## Elution Volume and Recovery Rate Comparing with Other Brands Products.

MonoSpin® ProA requires only 100 µL elution buffer to obtain a recovery rate of at least 90% IgG. On the other hand, other brands products requires 400 µL or more elution buffer with a recovery rate of 70% IgG.




## Removal of Preservatives in anti-body solutions



MonoSpin ProA/ProG enables you to remove proteins such as BSA and Gelatin in anti-body solutions without dilution.

# MonoSpin ProA, MonoSpin ProG

## Recovery Rate and Reproducibility of IgG from medium cultured CHO cells with MonoSpin® ProA 96 Well Plate



Sample volume : 150  $\mu$ L  
Elution volume : 150  $\mu$ L  
Recovery rate : 90 % (CV 3.1 %)  
IgG concentration : 1.3 mg/mL

## Purification of multiple antibodies using MonoSpin L and ProA

16 mg of antibody samples can be purified under vacuum or centrifugation complying with the following procedure.

### Procedure

1. Apply 5 mL of equilibration buffer.



2. Apply sample(Max. 8 mL) after filtration through 0.2.  $\mu$ L filtration.



3. Apply 5mL of washing buffer.



4. Apply 5mL of elution buffer.


Centrifugal force at each step: 1,500 x g, 2 min

\*MonoSpin ProA/G buffer kit was used.



# MonoSpin 96 well plate

MonoSpin 96WP is chosen for purification of multi-analytes. Immobilized silica monolith gel provides similar performance as MonoSpin series.



## Features

- Silica monolith gel are used for MonoSpin 96WP
- Performs perfectly for both Vacuum and centrifugation
- Rapid pretreatment of biological sample

## Application

- Purification and fractionation of peptides
- Recovery and purification of proteins
- Purification of sugar chains
- Purification of organic acid
- Recovery of drugs from biological samples (urine, serum, plasma)
- Purification of catecholamines

Description	Quantity	Cat.No.
MonoSpin 96WP C18	1	5010-21900
MonoSpin 96WP NH2	1	5010-21901
MonoSpin 96WP PBA	1	5010-21902
MonoSpin 96WP SAX	1	5010-21903
MonoSpin 96WP SCX	1	5010-21904
MonoSpin 96WP Amide	1	5010-21905
MonoSpin 96WP CBA	1	5010-21906
MonoSpin 96WP C18-CX	1	5010-21907
MonoSpin 96WP C18-AX	1	5010-21908



# MonoSpin 96 well plate

## 96 Deep Well Plate




- It is made of polypropylene and has excellent heat and cold tolerance as well as solvent tolerance.
- Superhydrophilic surface treatment can suppress nonspecific adsorption of proteins and peptides.
- Low adsorption(LB type) can prevent from adsorption of peptides and proteins.


Description	Material	Cat. No.	Qty	Cat.No.
MS Plate	polypropylene	SMST-0201	50	6045-00201
MS Plate Low adsorption <LB type>	polypropylene (hydrophobic polymer)	SMST-0801-LB	15	6045-00203

## GL Sticker for 96 well plate

### Evapo Less Slit



Enlarged view



- Sticker closes automatically after each apply.
- Adhesive-free on top of the sticker to prevent from contamination.
- Can be operated under  $-80^{\circ}\text{C} \sim 100^{\circ}\text{C}$

### Sealing Sticker



Enlarged view



- High durability against organic solvent
- High air leakage efficiency
- It can be used for storage of samples as it can bear up to  $-80^{\circ}\text{C}$ .

Description	Material	Cat. No.	Qty	Cat.No.
MS Plate	polypropylene	SMST-0201	50	6045-00201
MS Plate Low adsorption <LB type>	polypropylene (hydrophobic polymer)	SMST-0801-LB	15	6045-00203

# Ordering information

## MonoSpin type S

Description	Qty	Cat.No.
MonoSpin C18	50	5010-21700
	100	5010-21701
MonoSpin C18 FF	50	5010-21670
	100	5010-21671
MonoSpin Ph	50	5010-21733
	100	5010-21734
MonoSpin C18-AX	50	5010-21735
	100	5010-21736
MonoSpin C18-CX	50	5010-21731
	100	5010-21732
MonoSpin SAX	50	5010-21720
	100	5010-21721
MonoSpin SCX	50	5010-21725
	100	5010-21726
MonoSpin NH2	50	5010-21710
	100	5010-21711
MonoSpin CBA	50	5010-21729
	100	5010-21730
MonoSpin Amide	50	5010-21727
	100	5010-21728
MonoSpin PBA	50	5010-21715
	100	5010-21716
MonoSpin TiO	50	5010-21705
	100	5010-21706
MonoSpin Trypsin HP [ KEEP COOL ]	30	7510-11302
MonoSpin ME	50	5010-21737
	100	5010-21738
MonoSpin Phospholipid	50	5010-21698
	100	5010-21699



MonoSpin Type S



Recovery tube  
(1.7mL)



Liquid waste tube

## MonoSpin type S Trial kit

Trial kits and custom kits are shipped with various types of columns packaged for initial method development.

Description	Content	Cat.No.
MonoSpin Trial Kit 1	C18, TiO, SCX, SAX 10 each	5010-21740
MonoSpin Trial Kit 2	C18, Amide, CBA, NH2 10 each	5010-21741
MonoSpin Trial Kit 3	SCX, SAX, CBA, NH2 10 each	5010-21742

## MonoSpin type L

Description	Qty	Cat.No.
MonoSpin L C18	30	7510-11320
MonoSpin L SAX	30	7510-11321
MonoSpin L SCX	30	7510-11322
MonoSpin L NH2	30	7510-11323
MonoSpin L CBA	30	7510-11324
MonoSpin L ME	30	7510-11325
MonoSpin L Phospholipid	30	7510-11326

MonoSpin L type Product name



MonoSpin L type

## MonoSpin 96 well plate

Description	Qty	Cat.No.
MonoSpin 96WP C18	1	5010-21900
MonoSpin 96WP NH2	1	5010-21901
MonoSpin 96WP PBA	1	5010-21902
MonoSpin 96WP SAX	1	5010-21903
MonoSpin 96WP SCX	1	5010-21904
MonoSpin 96WP Amide	1	5010-21905
MonoSpin 96WP CBA	1	5010-21906
MonoSpin 96WP C18-CX	1	5010-21907
MonoSpin 96WP C18-AX	1	5010-21908

## MonoSpin ProA, MonoSpin ProG

Description	Qty	Cat.No.
MonoSpin ProA Column [ KEEP COOL]	10	7510-11310
MonoSpin ProG Column [ KEEP COOL]	10	7510-11311
MonoSpin ProA 96 Well plate [ KEEP COOL]	1	7510-11312
MonoSpin ProG 96 Well plate [ KEEP COOL]	1	7510-11313
MonoSpin L ProA [ KEEP COOL]	4	7510-11314
MonoSpin L ProG [ KEEP COOL]	4	7510-11315
MonoSpin ProA/G buffer kit [ KEEP COOL]	-	7510-11316

\*Various reagents required for purification is already attached

\*GL-SPE miniature suction manifolds are recommended for vacuum aspiration of 96-well plates

\*Centrifugal adaptor is attached to L type

## Global Solution

**GL Sciences**  
<https://www.gls.co.jp>

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

### **GL Sciences, Inc. Japan**

22-1 Nishishinjuku 6-Chome  
Shinjuku-ku, Tokyo,  
163-1130, Japan  
Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@gls.co.jp](mailto:world@gls.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

### **GL Sciences B.V.**

De Sleutel 9  
5652 AS Eindhoven  
The Netherlands  
Phone: +31 (0)40 254 95 31  
Email: [info@glscienc.es.eu](mailto:info@glscienc.es.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

### **GL Sciences (ShangHai) Ltd.**

Tower B, Room 2003,  
Far East International Plaza,  
NO,317 Xianxia Road,  
Changning District.  
Shanghai, China P.C. 200032  
Phone: +86 (0)21-6278-2272  
Email: [contact@glscienc.es.com.cn](mailto:contact@glscienc.es.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

### **GL Sciences, Inc. USA**

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
Phone: 310-265-4424  
Fax: 310-265-4425  
Email: [info@glscienc.esinc.com](mailto:info@glscienc.esinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

### **International Distributors**

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>

