



Column Selection Overview

Separation Mode Selection for Analyte Polarity



Column Selection for Analyte Properties



Water Soluble Vitamins	Phase	Detector
Vitamin C	NH ₂ , SIL, Amide	UV
Vitamin B ₁	ODS	UV, (FL)
Vitamin B ₂	ODS	UV, FL
Vitamin B ₆	ODS	UV
Niacin	ODS	UV
Folic acid	ODS	UV
Pantothenic acid	ODS	UV

Fat-Soluble Vitamin	Phase	Detector
Vitamin A	ODS	UV, FL
Vitamin D	ODS, NH ₂	UV
Vitamin E	ODS	UV, FL
Vitamin K	ODS	UV, (FL)

* (FL): Post-column derivatization

